



Fluorometer **LB-10FFM**

Index

Sr.no	Title	Page no
1.	Safety Measures	2
2.	Introduction	3
3.	Features	3
4.	Specifications	4
5.	Applications	4
6.	Instrument Introduction	5
7.	Operations	6
8.	Software Operations	9
9.	Maintenance	18
10.	Troubleshooting	18
11.	Accessories	18

1. Safety Measures

1.1 Operation Information

Users should know the working principle of the instrument before operation. Read this manual carefully.

1.2 Safety Instructions

The operation, maintenance, and repair of the instrument should comply with the basic guidelines and the remarked warning below.

If the users don't comply with them, it will affect the instrument.



Indoor used instrument.



Read the Manual carefully before operation. Only trained personnel can operate the instrument.



Turn off the power in case of stop working, unplug the power connector from the socket in case of long-term non-use, and cover the instrument with a cloth to avoid dust.



Unplug the power connector, contact your maintenance personnel who were trained in case of the following:

- There is some liquid flowing into the Instrument.
- The instrument gets wet through rain or water.
- Can't work normally, especially with odd sounds or smells.
- The instrument falls or comes out of the shell, damaged.

2. Introduction

Fluorometer LB-10FFM offers a wide concentration range down to 0.5pg/ μ L, ensuring high-sensitivity fluorescence detection. It features a dual-channel detection system for simultaneous measurement of two fluorescence signals. It is equipped with a monochrome LED light source and a high-sensitivity photodiode detector. Our fluorometer includes a 7-inch touchscreen display, enabling precise and reliable molecular analysis.

3. Features

- High-sensitivity detection
- Fast signal processing
- User-friendly interface
- Multiple calibration options
- Multiple excitation wavelengths
- Multi-format data export

4. Specifications

Model No.	LB-10FFM
Sample Capacity	1
Sample Volume	1 to 20 μ l
Tube Type	0.5ml PCR tube
Channel Number	2
Detection Time	3s
Repeatability	< 1.5%
Concentration Range	0.01ng/ μ l to 120ng/ μ l dsDNA HS 0.05ng/ μ l to 240ng/ μ l Oligo 0.2ng/ μ l to 2000ng/ μ l dsDNA BR 0.1mg/ml to 20mg/ml protein BR
Excitation Wavelength	470/625 (standard) 365/525 (optional)
Emission Wavelength	525/690 (standard) 460/620 (optional)
Calibration Method	2 or 3-point calibration
Response Range	Five orders of magnitude
Linearity	$R^2 > 0.995$
Light Source	Monochrome LED
Detector	Photodiode
Display	7-inch touchscreen display
Programs Stored	10000
Data Export	U Disk
Data Format	CSV, PDF
Data Interface	USB Type-A ($\times 2$), USB Type-B ($\times 1$)
Voltage	DC12V 2A
Power Consumption	4.5W
Power Supply	100-240V 50/60Hz
Dimensions (W \times D \times H)	296 \times 161 \times 55 mm
Packaging Dimension	390 \times 290 \times 190 mm
Net weight	1 kg
Gross Weight	2 kg

5. Applications

Fluorometer LB-10FFM is used for the qualitative and quantitative analysis of nucleic acids, proteins, and biomolecules in molecular biology and genetic studies. It is widely used in research institutes, pharmaceutical companies, clinical laboratories, and forensic science facilities.

6. Instrument Introduction

Instrument Structure:



Figure-1

Note: USB B for Program burn port, USB A for Printer, and USB disk (data export)

7. Operations

7.1 Working Conditions

- 1) **Ambient temperature:** 5°C – 35°C
- 2) **Relative humidity:** ≤70%
- 3) **Voltage:** DC12V 2A

7.2 Basic Operation



Figure-2 Insert the sample tube into the detection chamber

Insert the sample tube into the sample chamber and close the lid when setting up a standard curve or detecting.

Note:

- 1) Don't pinch tubes with standard or samples with fingers for a long time before detection, as body temperature will affect the solution in the tubes.
- 2) Do close the chamber lid to avoid interference with fluorescence from light during detection.
- 3) In case of tubes with scale lines, place tubes with the same direction for each measurement for the standard curve and detections to gain accurate and stable results.
- 4) The chamber temperature will heat up if it is closed for a long time, which will affect fluorescence sensitivity. When multiple detections are required from one PCR tube, keep the tube in a box free from light for 30 seconds before the next measurement.

7.3 Reagent preparation

Different preparation procedures for different kinds of reagents. Refer to the instructions that come with the reagent.




Figure-3 DNA/RNA/Protein detection reagent preparation

Reagent preparation steps:

- 1) **Standard preparation:** Get two 0.5ml PCR tubes, add each tube with 10µl of standard 1 and standard 2, and add 190µl of working solution respectively. Vortex mix both tubes for 2-3s, incubate at room temperature away from light for 2min or 15min. Then, start to set up the standard curve.
- 2) **Sample preparation:** Add the sample into a new tube with a volume 1-20µl, then add some working solution into the tube till the total volume to 200µl, vortex mix the tube for 2- 3s, incubate the tube at room temperature away from light for 2min or 15min, then start to detect. The steps above are only for DNA/RNA/Protein measurement; DNA/RNA requires 2 standard samples, while Protein requires 3.




Figure-4 Protein BR detection reagent preparation

Preparation steps are different between “Protein BR” and “DNA/RNA/protein”:
Standard preparation: Get two 0.5ml PCR tubes for S.1 and S.2 with 20µl standard 1 and standard 2, add 30µl reagent and 150µl buffer to each tube, respectively, to make a total volume of 200µl. Vortex mix each tube for 5s, incubate at room temperature for 10 min, then start to detect.

Note:

- 1) Room temperature of 22-25°C is recommended, according to measured data and relative data.
- 2) The DNA/RNA reagent fluorescence stability time is 3 hours, oligo reagent fluorescence stability time is 30 minutes at room temperature after incubation.

8. Software Operations

8.1 Home page













Figure-5

Icon	Description
	Calculator for reagent and concentration range
	Setting, for time settings and calibration
	DNA detection button
	RNA detection button
	Oligo detection button

Fluorometer LB-10FFM

 Protein	Protein detection button
 Fluorometer	This function is to detect sample fluorescence intensity. Fluorescence intensity is proportional to concentration, it's only used to detect relative concentration comparison for unknown samples.
Proj: <input type="text" value="YNXM"/>	Project name, tap to set for the project
ID: <input type="text" value="YNID"/>	ID name, tap to set
Type: <input type="text" value="Protein"/>	Tap to choose the sample type
Volume (μ l) <input type="text" value="10"/>	Tap to set sample volume in a 0.5ml PCR, can adjust the value using by buttons of + or -, can also input a value directly
Unit <input type="text" value="ng/<math>\mu</math>L"/>	Unit of concentration for different samples
Conc. (ng/ μ l) 4.89	Detection result
<input type="button" value="Detect"/>	Tap this button to start detection. Insert the sample tube into the chamber, close the lid, and tap this button to detect
 Curve	Tap this button to set up the standard curve
 Report	Tap this button to check the data

8.2 Set up a standard curve




Figure-6 Standard Curve



Figure-7 Set parameters for samples

Tap button “**Curve**” to set up the standard curve, utilize 2 or 3 standard samples with certain concentration to set up a standard curve before detection unknow samples.

Fluorometer LB-10FFM

Icon	Description															
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #0056b3; color: white;"> <th style="text-align: left;">Sample</th> <th style="text-align: left;">Conc.(ng/μl)</th> <th style="text-align: left;">FRU</th> </tr> </thead> <tbody> <tr> <td>Standard 1</td> <td>10</td> <td>80413</td> </tr> <tr> <td>Standard 2</td> <td>10</td> <td>80413</td> </tr> <tr> <td colspan="3"><hr/></td> </tr> <tr> <td colspan="3">Last calibration 2022-04-03 17:12:34</td> </tr> </tbody> </table>	Sample	Conc.(ng/μl)	FRU	Standard 1	10	80413	Standard 2	10	80413	<hr/>			Last calibration 2022-04-03 17:12:34			Showing the standard concentration, fluorescence intensity, and calibration time after setting the standard curve.
Sample	Conc.(ng/μl)	FRU														
Standard 1	10	80413														
Standard 2	10	80413														
<hr/>																
Last calibration 2022-04-03 17:12:34																
<div style="background-color: #0056b3; color: white; padding: 5px; text-align: center; width: fit-content; margin: auto;">Calibrate</div>	Tap to set the standard curve															
Conc. <input style="width: 100px;" type="text" value="10"/>	Input the standard sample concentration value															
Unit <input style="width: 100px;" type="text" value="ng/μL"/>	Tap to choose the Unit for a standard sample															
<div style="background-color: #0056b3; color: white; padding: 5px; text-align: center; width: fit-content; margin: auto;">Read standards</div>	After inputting the concentration value and Unit for standards, insert the standard sample into the chamber, close the lid, and tap this button to read the fluorescence value															

Different quantity of standard required for different type samples, after all standard samples detection, a standard curve shall be set up. Return to the home page to start detection.

8.3 Report








Figure-8 This page is for detection data

Fluorometer LB-10FFM

Tap “**Report**” on the home page.

Icon	Description						
Proj: <input type="text" value="YNXM"/> ID: <input type="text" value="YNID"/> Type: <input type="text" value="Protein"/>	Choose names for Proj, ID, Type to select data						
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #0070C0; color: white;"> <th style="text-align: left;">No.</th> <th style="text-align: left;">Detection time</th> <th style="text-align: left;">C(ng/μl)</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>2022-09-03 17:06:02</td> <td>7.92</td> </tr> </tbody> </table>	No.	Detection time	C(ng/μl)	1	2022-09-03 17:06:02	7.92	Display the detection time and result
No.	Detection time	C(ng/μl)					
1	2022-09-03 17:06:02	7.92					
<input type="button" value="First page"/> < <input type="text" value="5"/> > <input type="button" value="Last page"/>	Choose a different page						
	Insert U disk to export data						
	Tap to Delete data						
	Tap to print data after connecting a printer						
	Tap to enter calculator page						

8.4 Fluorometer




Figure-9 Fluorometer page

Fluorometer LB-10FFM

Icon	Description
Type: ● Blue470nm ▾	Tap to choose excitation light source, red, blue, or blue&red
FRU: 5432	Fluorescence value
Detect	Tap to detect the fluorescence value

8.5 Calculator




Figure-10 Calculator

This function is to calculate the volume of samples and standards.

Icon	Description
Samples: 10	Quantity of samples
Standards: 10	Quantity of samples
<input checked="" type="checkbox"/> Include overage (1 extra tube)	Option of the extra tube. Tick to add one extra one, untick to cancel
<ul style="list-style-type: none"> • Dye 0 μl <hr style="border: 0.5px dashed gray;"/> • Buffer 10 μl <hr style="border: 0.5px dashed gray;"/> • Total volume 10 μl 	Display volumes of dye, buffer, and total volume

Fluorometer LB-10FFM

<input type="button" value="Done"/>	Tap this button to start the calculation
-------------------------------------	--

Note: The Reagent calculator is not suitable for protein BR.

8.6 Data calculator

Molarity calculator

Desired units: →

Molecular weight: g/mol

Auto-populate DNA length

Samples	Original conc.(ng/μL)	Length (bp)	Mol.Conc.(μM)
S1	9999	<input type="text" value="10"/>	9999
S2	9999	<input type="text" value="10"/>	9999
S3	9999	<input type="text" value="10"/>	9999
S4	9999	<input type="text" value="10"/>	9999
S5	9999	<input type="text" value="10"/>	9999

Figure-11 Molarity calculator

Normalization

Calculator: Mol Conc. Weight

Final mol. concentration: μM









Final volume: μL

Samples	Original conc.(ng/μL)	Volume.(μL)	Diluent vol.(μL)
S1	9999	<input type="text" value="10"/>	9999
S2	9999	<input type="text" value="10"/>	9999
S3	9999	<input type="text" value="10"/>	9999
S4	9999	<input type="text" value="10"/>	9999
S5	9999	<input type="text" value="10"/>	9999

Figure-12 Normalization

Fluorometer LB-10FFM

On the data page, tap the “**Data**” button, which will bring you to the data calculator page, including the Molarity calculator and Normalization. The molarity calculator converts molecular weight to Molarity concentration. Normalization is used to convert several different concentration samples to the same molarity, concentration, or weight.

Icon	Description
	Set the unit before converting
	Set molecular weight
	Length (bp) is for the number of base pairs
	Tap to calculate the molarity concentration
	Choose the type of normalization
	Set concentration and volume for normalization
	Original concentration, volume required by normalization, and diluent volume.
	Tap to start normalization

8.7 Data and Time







Figure-13 Setting page

Tap the button in the upper right corner of the Home page to set Date and Time.

Icon	Description
	Input values directly to set Year, Month, Day
	Input values directly to set Hour, Minute, Second
	Tap this button to confirm settings

9. Maintenance

- 1) Use a soft, lint-free cloth dampened with clean water to clean the sample pedestal.
- 2) Do not use alcohol or solvent-based cleaners on the sample pedestal, as they may cause damage.
- 3) Dry the pedestal gently with a clean, dry cloth or allow it to air dry.
- 4) Wipe the instrument's exterior using a soft cloth and a mild cleansing cream or gentle detergent solution.
- 5) Avoid abrasive cleaners, scouring pads, or rough materials that can scratch the surface.
- 6) Prevent any liquid from entering the openings or ports of the instrument.
- 7) Keep cleaning materials clean and store them separately to avoid contamination.
- 8) Always turn off and unplug the instrument before cleaning.

10. Troubleshooting

Trouble	Cause	Solution
Unable to turn on	No electricity	Check the input port connection
Light source defect	Poor module connection	Check the instrument by shaking it gently to see if the module is fixed well or not
Unstable detection results	The PCR tube is not placed properly	Press the tube down again gently to a better position.
Inaccurate detection results	A standard curve is not good	Prepare standard samples and set up the standard curve again according to the steps strictly

11. Accessories

Standard Accessories

S. No	Accessory Name
1	U Disk



Labotronics Scientific. 1007 N Orange St., Wilmington, DE 19801, USA.
info@labotronics.com | www.labotronics.com