


Digital Colony Counter LB-20CC



Overview :

Digital Colony Counter LB-20CC is a tabletop counting unit with lamp power 16 W. Large diameter of counting plate ensures fast and reliable counting. Features three digit microbes counting capacity with LED light source for accurate counting.

Specifications :

Lamp power	16 W
Display digit	LED (0 to 999)
Dish diameter	? 115 mm
Light source	LED
Total power consumption	20 W
Power	AC 220 V,50 Hz
Dimensions	255 x 210 x 160 mm
Weight	4 kgs

Features :

1. Digital counting unit
2. LED light source for easy viewing of colonies
3. Counting plate diameter - ?115 mm
4. Counting capacity - upto three digit

Applications :

Used in healthcare, life science, molecular biology, pharmaceutical, clinical and research laboratories, microbiology, to estimate colonies of bacteria or other microorganism on agar, slide, gel or petri dish.



Labotronics Scientific. 1007 N Orange St., Wilmington, DE 19801, USA
info@labotronics.com | www.labotronics.com