



# Ducted Fume Hood LB-41DFH



## Overview :

Ducted Fume Hood LB-41DFH is designed to provide a safe and efficient working environment by effectively capturing and removing hazardous fumes and contaminants. Our fume hood operates within an airflow velocity range of 0.3 to 0.5 m/s Utilized with an aluminum alloy handle facilitates a seamless flow of air into the workspace, ensuring uninterrupted and effective airflow. Convenient handling with worry-free maintenance.

## Specifications :

Maximum Opening	780 mm
Airflow velocity	0.3~0.5 m/s
Illumination	>400LUX, LED lighting
Diameter of exhaust outlet	250 mm
Power supply	220/380V,50/60Hz, single or three phase.
Noise	