



Ducted Fume Hood LB-42DFH



Overview :

Ducted Fume Hood LB-42DFH is upgraded to effectively handle the removal of hazardous or volatile chemicals, ensuring a contamination-free environment. Operating with an airflow velocity ranging from 0.3 to 0.5 m/s maintains continuous airflow throughout the workstation. Equipped with an energy-saving microprocessor controller that controls humidity and temperature to provide accurate outcomes.

Specifications :

Maximum Opening	780 mm
Airflow velocity	0.3~0.5 m/s
Illumination	>400LUX, LED lighting
Diameter of exhaust outlet	250 mm
Power supply	220/380V,50/60Hz, single or three phase.
Noise	