



Ducted Fume Hood LB-43DFH



Overview :

Ducted Fume Hood LB-43DFH featuring an aluminum alloy handle, facilitates smooth airflow into the workspace, ensuring consistent and efficient ventilation. Engineered to create a safe working environment, it effectively captures and eliminates hazardous fumes and contaminants. Operating within an airflow velocity range of 0.3 to 0.5 m/s, it ensures optimal performance. With convenient handling and hassle-free maintenance, it offers ease of use for users.

Specifications :

Maximum Opening	780 mm
Airflow velocity	0.3~0.5 m/s
Illumination	>400LUX, LED lighting
Diameter of exhaust outlet	250 mm
Power supply	220/380V,50/60Hz, single or three phase.
Noise	