



Ductless Fume Hood LB-20DLFH



Overview :

Ductless Fume Hood LB-20DLFH is a benchtop unit comes with safe and efficient air filtration and exhaust system. Features low noise design with high quality DC centrifugal blower to achieve quiet working environment. With 3-sided acrylic transparent observation window for multi-angle observation and LCD touch screen display, ensure convenient operation.

Specifications :

Air velocity range	0.4 to 0.6 m/s
Maximum opening	650 mm
Operation surface height	690 mm
Blower	DC centrifugal blower
Filter	HEPA filter × 1 pc
Lamp	LED lamp
Display	7-inch color touch screen
Operation chamber dimension (W×D×H)	780×540×775 mm
Dimension (W×D×H)	800×620×1130 mm
Packaging dimension (W×D×H)	980×980×720 mm
Gross weight	80 kg

Features :

1. A benchtop ductless fume hood ideal to use for different applications
2. Safe and efficient air filtration system and exhaust system, reliable operation
3. Chemical and HEPA filter with fan system for effective sterile environment
4. Low noise design with high quality DC centrifugal blower for high efficiency
5. 3-sided acrylic transparent observation window for multi-angle observation
6. LCD touch screen display for real-time display of different operation parameters
7. Not require duct system installation, allow flexible placement at any location

Applications :

Ductless fume hood is used to remove hazardous fumes, vapors and particles without any environmental disturbance across laboratories and research field where outside ventilation cannot be achieved.

Optional Accessories:

- Base stand
- Base cabinet
- Active carbon filter
- LED lamp