



Narcotic Chemical Storage Cabinet LB-15NSC



Overview :

Narcotic Chemical Storage Cabinets LB-15NSC are corrosion and a fire resistant units. With a double locking system the cabinets are designed to store and manage narcotics and medicines that require strict temperature control and high level protection.

Specifications :

Capacity (litres)	340
Loading capacity	100 kgs
Number of Shelves	2
Number of doors	2
Door style	Manual
Dimensions (H × W × D)	1650 × 1090 × 860 mm
Colour	Grey
Weight	190 kgs

Features :

1. Double locking system (Digital lock + Key lock)
2. Durable lead free epoxy powder coated material for chemical resistance
3. Dual ventilation holes with flame arresters
4. Adjustable legs and grounding wires

Applications :

Used in chemical industries, laboratories, pharmacies, hospitals, scientific research institutions, etc.

Optional Accessories:

1.	PP Shelves
----	------------