



Paraffin Dispenser LB-10PD



Overview :

Paraffin dispenser LB-10PD is a microprocessor controlled unit for liquid paraffin dispensing. It is characterized with reservoir for homogeneous heat dissipation. The modulated feature ensures retention of molecular characteristics of paraffin.

Specifications :

Paraffin reservoir capacity	5 L
Temperature range	30 °C ~ 90 °C
Temperature increments	1 °C
Heating plate dimensions	310 mm × 150 mm
Power supply	AC 220 V; 50 Hz AC 110 V; 60 Hz
Dimensions	310 × 350 × 410 mm
Weight	13 kgs

Features :

1. Digital LED screen for visual monitoring of set parameters
2. Thermal shock prevention by homogeneous heat dissipation
3. Adjustable temperature setting for quick paraffin melting
4. Tap assembly for proper paraffin pour
5. Anti-corrosive structural unit

Applications :

Used for quick access to liquid paraffin inn histopathology labs, pathological labs, zoological and botanical labs, medical and veterinary institutes, scientific research and development, etc.