



# Table Top Autoclave LB-10TBA



## Overview :

Class B Table top autoclave LB-10TBA is a microcomputer control system with three times better pre-vacuum drying function. The built-in steam generator provides rapid cycles of saturated steam. It is characterized with automatic cooling air and steam exhaustion function at the end of the sterilizing process.

## Specifications :

Capacity	8 L
Sterilization room size	Φ 170 × 360 mm
External dimensions	410 × 530 × 400 mm
Working temperature	134 °C
Maximum pressure	0.22 MPa
Power	1500 W
Voltage	110 V; 60 Hz 220 V; 50 Hz
Weight	42 kgs

## Features :

1. Malfunction detection (over temperature, over voltage and low water level) by indication
2. European Class B standard with protection function over temperature, pressure, low water level and anti-burning
3. Buzzing alarm indication and automatic process end on completion of sterilization
4. Interlocking system to prevent spillage of water from top tank
5. High quality stainless steel liner with sealed cover
6. Equipped with standard test interface
7. Digital display with touch screen panel for ease of operation
8. Automatic door switch and consists quality test
9. Provides three times per-vacuum and saturated steam through built in generator

## Applications :

Used in sterilizing dental tools, veterinary medicines, biological products, utensils, medical appliances, culture media, etc. in the fields of dentistry, biology, biochemistry, forensic, pharmaceuticals, cytology, immunology, genetics, bioengineering, health and sanitary, etc.